

---

**Product Data Sheet**


---

Product Name: NRIP3 Human  
 Cat. No.: GP24065  
 Batch No.: 1

**Product Data**

Purity	>98%	Source	E.coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	C11orf14; NY-SAR-105; Nuclear receptor-interacting protein 3; Sarcoma antigen NY-SAR-105; NRIP3.		
Amino Acid Sequence	MGSSHHHHHH SGLVPRGSH MGSMFYSGLL TEGGRKETDM REAASLRQQR RMKQAVQFIH KDSADLLPLD GLKKGSSKD MQPHNILQRR LMETNLSKLR SGPRVPWASK TNKLNQAKSE GLKKSEEDDM ILVSCQCAGK DVKALVDTGC LYNLISLACV DRLGLKEHVK SHKHEGEKLS LPRHLKVVGQ IEHLVITLGS LRLDCPAAVV DDNEKNLSLG LQTLRSLKCI INLDKHRLIM GKTDKEEIPF VETVSLNEDN TSEA.		
Formulation	NRIP3 protein solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 20% glycerol and 1mM DTT.		

**Introduction**

Nuclear Receptor-Interacting Protein (NRIP3) is a 241 amino acid protein which is encoded by a gene that maps to human chromosome 11. Diseases associated with NRIP3 include sarcoma, and acute myeloid leukemia.

**Stability**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Background**

NRIP3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 264 amino acids (1-241 a.a.) and having a molecular mass of 29.4kDa. NRIP3 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

---

**Product Data Sheet**

---

***Caution: Product has not been fully validated for medical applications. For research use only.***

***Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com***

***Address: 10292 Central Ave. #205, Montclair, CA, USA***