

---

**Product Data Sheet**

---

Product Name: Nucleosome Bovine

Cat. No.: GP24080

Batch No.: 1

**Product Data**

Purity	>98%	Source	Bovine tissues.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	Nucleosome Bovine is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.		

**Introduction**

Nucleosome core proteins are a group of proteins which are necessary for DNA packaging in eukaryotes. Every core contains an octamer with 2 copies of each histone 2A, 2B, 3 and 4 around which an approx. 140 bp section of DNA is wound twice.

**Stability**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

**Background**

Nucleosome Bovine is purified from bovine tissues by proprietary protein-chemical techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA