

---

**Product Data Sheet**


---

Product Name: S100A8 Human, His-Myc  
 Cat. No.: GP24474  
 Batch No.: 1

**Product Data**

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	S100 Calcium Binding Protein A8; S100 Calcium-Binding Protein A8 (Calgranulin A); Migration Inhibitory Factor-Related Protein 8; Leukocyte L1 Complex Light Chain ; Urinary Stone Protein Band A; Calprotectin L1L Subunit; Cystic Fibrosis Antigen; Calgranulin-A; MRP8; CAGA; CFAG; P8; S100 Calcium Binding Protein A8 (Calgranulin A); S100 Calcium-Binding Protein A8; Calgranulin A; 60B8AG; CP-10; MA38; MRP-8; CGLA; L1Ag; MIF; NIF.		
Amino Acid Sequence	MGSSHHHHHSSGLVPRGSH MGSFMLTELEKALNSIIDVY HKYSLIKGNF HAVYRDDLKK LLETECPQYI RKKGADVWFK ELDINTDGAV NFQEFLILVIKMGVAAHKKS HEESHKEEQK LISEEDL.		
Formulation	The S100A8 solution (1mg/ml) contains 20mM Tris-HCl (pH8.0), 10% glycerol and 0.1M NaCl.		

**Introduction**

S100A8 is a part of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a broad range of cells, and participate in the regulation of cellular processes such as cell cycle progression and differentiation. S100A8 plays a role in the inhibition of casein kinase and as a cytokine. S100A8 altered expression is related with cystic fibrosis disease. S100A8 is a calcium-binding protein that has antimicrobial activity against bacteria and fungi. S100A8 is crucial for resistance towards invasion by pathogenic bacteria. S100A8 up-regulates transcription of genes that are under the control of NF-kappa-B. S100A8 plays a role in the development of endotoxic shock in response to bacterial lipopolysaccharide. S100A8 endorses tubulin polymerization and promotes phagocyte migration and infiltration of granulocytes at sites of wounding. S100A8 takes part as a

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---

## Product Data Sheet

---

pro-inflammatory mediator in acute and chronic inflammation and up-regulates the release of IL8 and cell-surface expression of ICAM1.

### Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Background

S100A8 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 127 amino acids (1-93) and having a molecular mass of 14.5kDa. S100A8 is fused to a 24 aa His-tag at N-terminus and to a 10 aa Myc-tag at C-terminus.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA