
Product Data Sheet

Product Name: SHFM1 Human
 Cat. No.: GP24571
 Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	SHFM1; Split Hand/Foot Malformation (Ectrodactyly) Type 1; DSS1; SHFD1; Deleted In Split Hand/Split Foot Protein 1; Split Hand/Foot Deleted Protein 1; Split Hand/Foot Malformation Type 1 Protein; Deleted In Split-Hand/Foot 1; 26S Proteasome Complex Subunit DSS1; ECD; SEM1; SHSF1; Shfdg1; Deleted In Split-Hand/Split-Foot 1; SHFDG1.		
Amino Acid Sequence	MGSSHHHHHH SGLVPRGSH MGSMSEKKQP VDLGLLEEDD EFEEFPAEDW AGLDEDEDAH VWEDNWDDDN VEDDFSNQLR AELEKHGYKM ETS.		
Formulation	SHFM1 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 10% glycerol and 0.1M NaCl.		

Introduction

26S proteasome complex subunit DSS1 (SHFM1) has been suggested as a candidate gene for the autosomal dominant form of the heterogeneous limb developmental disorder split hand/split foot malformation type 1. SHFM1 is a subunit of the 26S proteasome which plays a part in ubiquitin-dependent proteolysis. SHFM1 binds and stabilizes BRCA2 and is therefore involved in the control of R-loop-associated DNA damage and thus transcription-associated genomic instability. Furthermore, SHFM1 may have a role in the completion of the cell cycle. SHFM1 is a component of the TREX-2 complex (transcription and export complex 2), comprised of at least ENY2, GANP, PCID2, DSS1, and either centrin CETN2 or CETN3.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Background

SHFM1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 93 amino acids (1-70 a.a) and having a molecular mass of 10.7kDa (Molecular size on SDS-PAGE will appear higher).SHFM1 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA