
Product Data Sheet

Product Name: SNRPD Human
Cat. No.: GP24641
Batch No.: 1

Product Data

Purity	>98%	Source	Sf9 insect cells.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	SNRPD is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.		

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

Background

SNRPD Human Recombinant produced in SF9 insect cells, is a mixed protein complex containing the SNRPD1 (14,306 Dalton), SNRPD2 (14,773 Dalton) and SNRPD3 (14,739 Dalton) proteins. Each protein (SNRPD1+ SNRPD2+ SNRPD3) is individually fused to a -6xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA