

---

## Product Data Sheet

---

Product Name: RNP/Sm Bovine  
Cat. No.: GP24643  
Batch No.: 1

### Product Data

Purity	>98%	Source	Bovine tissues.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	RNP/Sm Bovine is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.		

### Introduction

U-snRNPs which are a pivotal part of spliceosomes, are created between small non-coding nuclear RNAs and a protein moiety. Whereas the 7 Sm proteins are conserved, the RNP antigens are exclusive for each U-snRNP subclass.

### Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

### Background

RNP/Sm Bovine is purified from bovine tissues by proprietary protein-chemical techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA