
Product Data Sheet

Product Name: Ubiquitin K48R Human
Cat. No.: GP24939
Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Ubiquitin; Ribosomal Protein S27a; CEP80; UBA80; UBCEP1; UBCEP80; HUBCEP80; RPS27A; Ubiquitin K48R.		
Formulation	Diluted in PBS plus 5% glycerol.		

Introduction

Ubiquitin is a highly conserved 76 amino acid protein expressed in all eukariotes. Ub is found either in free form or conjugated to proteins as monomer or as chain of ubiquitin molecules. The most well characterized consequence of polyubiquitination is substrate degradation, while mono-ubiquitinated proteins are not degraded. Ubiquitin attachment to protein substrate is a complex process involving a ubiquitin activating enzyme (E1), a ubiquitin conjugating enzyme (E2) and a ubiquitin protein ligase (E3). The first ubiquitin moiety is transferred to the e-NH₂ group of a Lys residue of the protein substrate to generate an isopeptide bond. In successive reactions, a poly ubiquitin chain is synthesized by processive transfer of additional activated moieties to Lys48 of the previously conjugated ubiquitin molecule. Ubiquitin K48R prevents the formation of poly ubiquitin chains via Lys48 linkages with mono ubiquitin molecules, avoiding the degradation of protein substrates.

Stability

Store vial at -20°C to -80°C . When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Background

Recombinant human ubiquitin featuring a Lys 48 to Arg48 mutation is useful for the reduction of poly-Ub chain length and conjugation rates. Ubiquitin K48R is expressed in E.coli and purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA