

---

**Product Data Sheet**

---

Product Name: Phl p 2  
Cat. No.: GP25142  
Batch No.: 1

**Product Data**

Purity	>98%	Source	Sf9 insect cells.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Pollen allergen Phl p 2; Allergen Phl p II; Phl p 2; PHLPII.		
Formulation	Phl p 2 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.		

**Introduction**

Pollen allergen Phl p2 (Phl p 2) causes an allergic reaction in humans.

**Stability**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.

**Background**

Recombinant Pollen allergen Phl p 2 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 14,734 Dalton. Phl p 2 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA