
Product Data Sheet

Product Name: Mumps
Cat. No.: GP25160
Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	Phosphate buffer, pH 7.4.		

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

Full length mumps viral nucleoprotein was expressed from E. coli, the recombinant protein was purified by affinity chromatography technique. The expressed recombinant mumps viral nucleoprotein contains 553 amino acids and a 6 x histidine tag was fused at its N-terminus, having a total Mw of 62kDa.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA