

---

## Product Data Sheet

---

Product Name: SARS Core (1-49)

Cat. No.: GP25167

Batch No.: 1

### Product Data

Purity	>98%	Source	
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	50mM Tris-HCl, 60mM NaCl and 50% glycerol.		

### Introduction

SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies.

### Stability

SARS Core protein is shipped at ambient temperature. Upon arrival, store at -20C.

### Background

The E.Coli derived 32 kDa recombinant protein contains the Nucleocapsid core protein 1-49 amino acids immunodominant regions.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA