
Product Data Sheet

Product Name: HCV NS5 Genotype-3

Cat. No.: GP25296

Batch No.: 1

Product Data

Purity	>98%	Source	
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	25mM Tris-HCl, 1mM EDTA, 1.5M urea and 50% glycerol.		

Introduction

HCV is a small 50nm, enveloped, single-stranded, positive sense RNA virus in the family Flaviviridae. HCV has a high rate of replication with approximately one trillion particles produced each day in an infected individual. Due to lack of proofreading by the HCV RNA polymerase, the HCV has an exceptionally high mutation rate, a factor that may help it elude the host's immune response. Hepatitis C virus is classified into six genotypes(1-6) with several subtypes within each genotype. The preponderance and distribution of HCV genotypes varies globally. Genotype is clinically important in determining potential response to interferon-based therapy and the required duration of such therapy. Genotypes 1 and 4 are less responsive to interferon-based treatment than are the other genotypes (2, 3, 5 and 6).

Stability

HCV NS5 Genotype-3 although stable at 4°C for 1 week, should be stored below -18°C .Please prevent freeze thaw cycles.

Background

The E.coli derived recombinant protein fused to a GST tag contains the HCV NS5 Genotype 3 immunodominant regions.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA