
Product Data Sheet

Product Name: HCV NS5 Genotype-4

Cat. No.: GP25299

Batch No.: 1

Product Data

| | | | |
|---------------------|---|--------------------|-------------------------|
| Purity | >98% | Source | |
| Physical Appearance | solid | Shipping Condition | Shipped with Ice Packs. |
| Formulation | 1.5M urea, 25mM Tris-HCl, pH-8, 0.2% Triton-X & 50% Glycerol. | | |

Introduction

HCV is a small 50nm, enveloped, single-stranded, positive sense RNA virus in the family Flaviviridae. HCV has a high rate of replication with approximately one trillion particles produced each day in an infected individual. Due to lack of proofreading by the HCV RNA polymerase, the HCV has an exceptionally high mutation rate, a factor that may help it elude the host's immune response. Hepatitis C virus is classified into six genotypes(1-6) with several subtypes within each genotype. The preponderance and distribution of HCV genotypes varies globally. Genotype is clinically important in determining potential response to interferon-based therapy and the required duration of such therapy. Genotypes 1 and 4 are less responsive to interferon-based treatment than are the other genotypes (2, 3, 5 and 6).

Stability

HCV NS5 Genotype-4 although stable at 4°C for 1 week, should be stored below -18°C .Please prevent freeze thaw cycles.

Background

The E.coli derived recombinant protein contains the HCV NS5 Genotype 4 immunodominant regions, amino acids 2212-2313. The protein is fused with GST tag at N-terminus.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA