
Product Data Sheet

Product Name: HBsAg adr CHO

Cat. No.: GP25346

Batch No.: 1

Product Data

Purity	Greater than 95.0% as determined by RP-HPLC.	Source	CHO.
Physical Appearance	Sterile Filtered clear solution.	Shipping Condition	Shipped with Ice Packs.
Formulation	Sterile Filtered solution containing 20mM PB and 154mM NaCl.		

Introduction

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV) S-gene. The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

Stability

HBsAg Should be stored at 4°C , DO NOT FREEZE.

Background

Recombinant HBsAg adr produced in CHO cells, the molecular weight is approximately 23-27kDa.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA