
Product Data Sheet

Product Name: Influenza-A H1N1 New Caledonia
Cat. No.: GP25488
Batch No.: 1

Product Data

Purity	>98%	Source	Baculovirus Insect Cells.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	The Recombinant H1N1 A/New Caledonia/20/99 solution 10mM Sodium phosphate, pH 7.2, 150mM NaCl.		

Introduction

H1N1 is a subtype specie of Influenza A virus. H1N1 Influenza Virus has mutated into various strains such as the Spanish Flu strain, mild human flu strains, endemic pig strains, and various strains found in birds. The Influenza A Virus is a globular particle about 100nm in diameter, sheathed in a lipid bilayer derived from the plasma membrane of its host. Studded in the lipid bilayer are two integral membrane proteins some 500 molecules of hemagglutinin ("H") and some 100 molecules of neuraminidase ("N"). Within the lipid bilayer are 3000 molecules of matrix protein and 8 pieces of RNA. Each of the 8 RNA molecules is associated with many copies of a nucleoprotein, several molecules of the three subunits of its RNA polymerase some "non-structural" protein molecules of uncertain function.

Stability

H1N1 A/New Caledonia/20/99 Recombinant should be stored at 4°C . Do not freeze!

Background

Recombinant Full-Length H1N1 New Caledonia/20/99 is glycosylated with N-linked sugars, produced using baculovirus vectors in insect cells and its Mw is approximately 72 kDa. The insect cells infected with A9440.1a recombinant baculovirus expressing recombinant H1N1 A/New Caledonia/20/99. H1N1 New Caledonia shows 90% similarity to the A/PR/8/34 amino acid sequence. The accession number is DQ508857.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA