
Product Data Sheet

Product Name: H3N2 Canine
 Cat. No.: GP25496
 Batch No.: 1

Product Data

Purity	>98%	Source	Sf9, Baculovirus cells.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Synonyms	Hemagglutinin		
Amino Acid Sequence	ADPNLPGNEN NAATLCLGHH AVPNGTIVKT ITDDQIEVTN ATELVQNSST GKICNNPHKI LDGRDCTLID ALLGDPHCDV FQNETWDLFV ERSNAFSNCY PYDVPDYASL RSIVASSGTL EFITEGFTWA GVTQNGGSGA CKRGPANGFF SRLNWLTKSG NTYPVLNVTM PNNNNFDKLY IWGVHHPSTN QEQTSLYIQA SGRVTVSTRR SQQTIIPNIG SRPLVRGQSG RISVYWTIVK PGDVLVINSN GNLIAPRGYF KMRIGKSSIM RSDAPIDTCI SECITPNGSI PNEKPFQNVN KITYGACPKY VKQNTLKLAT GMRNVPEKQT HHHHHH.		
Formulation	Canine H3N2 protein solution (0.5mg/ml) contains 10% glycerol & Phosphate Buffered Saline (pH 7.4).		

Introduction

H3N2 is a subtype of the influenza A virus. Its name derives from the forms of the two kinds of proteins on the surface of its coat, hemagglutinin (H) and neuraminidase (N). H3N2 exchanges genes for internal proteins with other influenza subtypes. H3N2 has tended to dominate in prevalence over H1N1, H1N2, and influenza B. H3N2 strain descended from H2N2 by antigenic shift, in which genes from multiple subtypes reassorted to form a new virus. Both the H2N2 and H3N2 strains contained genes from avian influenza viruses. H3N2 viruses are able to infect mammals and birds. In pigs, humans, and birds, the virus has mutated into many strains. Hemagglutinin(HA) binds to sialic acid-containing receptors on the cell surface, generating the attachment of the virus particle to the cell. HA has a vital part in the determination of host range restriction and virulence and is in charge of the diffusion of the virus into the cell cytoplasm by facilitating the fusion of the membrane of the endocytosis.

Stability

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Background

Canine H3N2 produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing 336 amino acids (18-344a.a.) and having a molecular mass of 36.9kDa. (Molecular size on SDS-PAGE will appear at approximately 40-57kDa). H3N2 is expressed with a 6 amino acid His tag at C-Terminus and purified by proprietary chromatographic techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA