
Product Data Sheet

Product Name: Dengue Envelope-1 45kDa

Cat. No.: GP25533

Batch No.: 1

Product Data

| | | | |
|---------------------|---|--------------------|-------------------------|
| Purity | >98% | Source | E.coli. |
| Physical Appearance | solid | Shipping Condition | Shipped with Ice Packs. |
| Formulation | Phosphate buffered saline with 25mM arginine. | | |

Introduction

Dengue fever is affected by 1 of 4 closely linked virus serotypes of the genus Flavivirus, family Flaviviridae. One might have dengue fever infected by the different serotype virus after the primary infection. Detection of particular antibodies to dengue viruses is used for the diagnosis in clinic. Recently, lateral flow rapid test products have become a most suitable and known method in the clinical diagnosis. Though, there is difficulty for manufacturers to have the dengue antigens with a whole coverage for Dengue IgG & IgM recognition for all 4 serotype infections as well as with colloid gold binding ability. 8 dengue antigens have been established for the lateral flow products which are well defined for their coverage of dengue IgG & IgM recognition. The researcher can select the products below based on their own specific application.

Stability

Dengue Envelope-1 45kDa although stable at 4°C for 1 week, should be stored below -18°C .Please prevent freeze thaw cycles.

Background

Recombinant Dengue Virus Subtype 1 Envelope 45kDa (43-413a.a.) produced in E.coli, contains important epitopes for dengue antibody recognition. This protein is fused to a 6xHis tag.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA