
Product Data Sheet

Product Name: Dengue-2 C-terminal

Cat. No.: GP25558

Batch No.: 1

Product Data

Purity >98%

Source E.coli.

Physical Appearance solid

Shipping Condition Shipped with Ice Packs.

Formulation PBS.

Introduction

Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

Stability

Dengue 2 Envelope Recombinant although stable at 4°C for 1 week, should be stored below -18°C .

Background

The Recombinant Dengue 2 Envelope Hydrophobic Domain e.coli derived contains 225 amino acids (228-485) from dengue 2 envelope. The protein is fused to a His tag at C-terminal and purified by standard chromatography techniques.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA