

---

## Product Data Sheet

---

Product Name: SARS Nucleocapsid, 1-49 a.a.

Cat. No.: GP26378

Batch No.: 1

### Product Data

Purity >98% Source Escherichia Coli.

Physical Appearance solid Shipping Condition withIcePacks

Formulation SARS Nucleocapsidprotein solution is supplied in PBS.

### Introduction

The SARS Nucleocapsid is a protein that has a crucial structural importance, by making the virions associates with the RNA genome and therefore create a nucleocapsid that has the specific helical form. There has been found 3 different domain structures to the nucleocapsid antigen. This protein contains 3 outer structural proteins, the membrane (M), envelope (E), and spike (S) proteins.

### Stability

Protein is shipped on ice packs. Upon arrival, Store at -20°C.

### Background

The E.coli derived recombinant protein contains the Nucleocapsid protein 1-49 amino acids immunodominant regions.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA