
Product Data Sheet

Product Name: SARS Membrane
Cat. No.: GP26380
Batch No.: 1

Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	withIcePacks
Amino Acid Sequence	S part a.a. 20 to 210		
Formulation	SARS Membraneprotein solution is supplied in PBS.		

Introduction

SARS Membrane protein is a membranal protein which plays a crucial role in the virus' life cycle. This protein has 3 outer cell structural antigens, the envelope (E), membrane (M) and spike (S) antigens. SARS Membrane protein is essential to the process of assembly, budding, creation of the envelope and pathogenesis.

Stability

Protein is shipped on ice packs. Upon arrival, Store at -20°C.

Background

The E.coli derived recombinant protein contains the Membrane/Matrix protein immunodominant regions fused to 6xHis tag at C-terminal.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA