

---

## Product Data Sheet

---

Product Name: SARS Spike (1-53, 90-115, 171-205)

Cat. No.: GP26381

Batch No.: 1

### Product Data

Purity >98% Source Escherichia Coli.

Physical Appearance solid Shipping Condition withIcePacks

Formulation SARS Spikeprotein solution is supplied in PBS.

### Introduction

SARS Associated Coronavirus Spike is a causative agent of SARS virus. Diseases caused by Coronaviruses are usually controlled by CD8 T cells. The spike antigen is one of the main surface proteins on Coronavirus, therefore one of the main proteins that are candidates for vaccines.

### Stability

Protein is shipped on ice packs. Upon arrival, Store at -20°C.

### Background

The E.coli derived recombinant protein contains the Spike (1-53, 90-115, 171-205 a.a.) protein immunodominant regions fused to 6xHis tag at C-terminal.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA