
Product Data Sheet

Product Name: SARS Spike (408-470, 540-573)

Cat. No.: GP26382

Batch No.: 1

Product Data

Purity >98% Source Escherichia Coli.

Physical Appearance solid Shipping Condition withIcePacks

Formulation SARS Spikeprotein solution is supplied in PBS.

Introduction

Severe acute respiratory syndrome (SARS) has a new causative agent identified as Coronavirus. SARS Associated Coronavirus Spike is one of the main surface antigens of Coronavirus, making it one of the main candidates for vaccines. Diseases caused by Coronaviruses are usually controlled by CD8 T cells.

Stability

Protein is shipped on ice packs. Upon arrival, Store at -20°C.

Background

The E.coli derived recombinant protein contains the Spike (408-470, 540-573 a.a.) protein immunodominant regions fused to 6xHis tag at C-terminal.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA