
Product Data Sheet

Product Name: SARS Spike Mosaic

Cat. No.: GP26383

Batch No.: 1

Product Data

Purity >98% Source Escherichia Coli.

Physical Appearance solid Shipping Condition withIcePacks

Formulation SARS Spike Mosaicprotein solution is supplied in1xPBS, pH 7.8.

Introduction

SARS or severe acute respiratory syndrome is a lethal kind of pneumonia generated by coronavirus. Coronaviruses are a family of viruses that causes a mild respiratory tract infection in a large type of creatures such as humans, pigs, cows, mice, cats etc. up to this day, the origin of the coronavirus is still unknown.

Stability

Protein is shipped on ice packs. Upon arrival, Store at -20°C.

Background

The E.coli derived recombinant protein contains the Spike Mosaic protein immunodominant regions 20-210 a.a. fused to 6xHis tag at C-terminal.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA