

---

## Product Data Sheet

---

Product Name: Yellow Fever Virus

Cat. No.: GP26449

Batch No.: 1

### Product Data

Purity	>98%	Source	Escherichia Coli.
Physical Appearance	solid	Shipping Condition	withIcePacks
Formulation	Recombinant Yellow Fever Virus protein solution contains PBS, 25mM K2CO3.		

### Introduction

Yellow fever virus is a part of Flaviviridaeis virus family. Yellow Fever Virus is a 40–50nm enveloped RNA virus, transmitted by mosquito. Recombinant yellow fever virus protein is derived from its envelope domain III.

### Stability

Upon arrival, Store at -20°C. Avoid multiple freeze-thaw cycles.

### Background

Recombinant Yellow Fever Virus produced in E. coli migrates at 12kDa. Recombinant Yellow Fever Virus is fused to a 6xHis tag at its C-terminus and purified by proprietary chromatographic technique.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA